

Client: Example Client ABC123 123 Test Drive Salt Lake City, UT 84108

UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB 4/10/1993

Sex:

Patient Identifiers: 01234567890ABCD, 012345

Male

Visit Number (FIN): 01234567890ABCD **Collection Date:** 01/01/2017 12:34

RBC Band 3 Protein Reduction in Hereditary Spherocytosis

ARUP test code 2008460

RBC Band 3 Protein Reduction in HS

Borderline

The findings are borderline. The possibility of hereditary spherocytosis cannot be ruled out completely. Recent transfusion (within the last 3 months) may affect results. Please correlate with clinical and other laboratory findings.

INTERPRETIVE INFORMATION: RBC Band 3 Protein Reduction in HS

This test can be used to confirm a suspected diagnosis of Hereditary Spherocytosis (HS). HS is a common inherited hemolytic anemia characterized by the presence of spherical erythrocytes (spherocytes). HS is diagnosed based on family history and clinical features, along with clinical laboratory tests, including peripheral smear examination, osmotic fragility (OF), flow cytometry, or by genetic testing (Hereditary Hemolytic Anemia Panel Sequencing. ARUP test code 2012052).

Band 3 (or solute carrier family 4 member 1, SLC4A1) is the most abundant transmembrane protein found in human red blood cells (RBC). Eosin-5-maleimide (EMA) dye binds to band 3 on intact RBC's. A reduction of fluorescence intensity will be seen in hereditary spherocytosis. This test by flow cytometry has been reported to have a sensitivity of 93 percent for a diagnosis of HS. Congenital Dyserythropoietic Anemia Type II, Southeast Asian Ovalocytosis and Hereditary Pyropoikilocytosis are rare disorders that may also show a positive result.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
RBC Band 3 Protein Reduction in HS	20-133-400883	5/11/2020 11:30:00 AM	5/13/2020 7:08:16 PM	5/15/2020 1:51:00 PM

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: