

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 2/19/1941
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Heat Shock Protein 70, IgG by Immunoblot

ARUP test code 2013590

Heat Shock Protein 70, IgG by Immunoblot

Positive * (Ref Interval: Negative)

INTERPRETIVE INFORMATION: Heat Shock Protein 70, IgG by Immunoblot

The presence of HSP70 IgG antibodies may be useful in predicting corticosteroid responsiveness in a subset of patients with autoimmune inner ear disease (AIED) characterized by idiopathic rapidly progressive sensorineural hearing loss (SNHL). HSP70 IgG antibodies are also associated with a number of autoimmune diseases and have also been reported in apparently healthy individuals. A negative result does not rule out response to treatment or to a diagnosis of AIED.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement D: aruplab.com/CS

VERIFIED/REPORTED DATES

| Procedure | Accession | Collected | Received | Verified/Reported |
|--|---------------|------------------|------------------|-------------------|
| Heat Shock Protein 70, IgG by Immunoblot | 19-031-140741 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical