

Patient: [REDACTED]  
 DOB: [REDACTED] Age: 30  
 Patient Identifiers: [REDACTED]  
 Visit Number (FIN): [REDACTED]

Sex: F

Client: ARUP Physician Services  
 321 TESTING ANSR EXTRACT  
 Salt Lake City, NY 84108

Physician: [REDACTED]


ARUP Test Code: 3002254

Collection Date: 09/19/2022

Received in lab: 09/19/2022

Completion Date: 09/20/2022


**Whole Peanut Specific IgE**



**REFERENCE INTERVAL: Allergen, Interpretation**

Less than 0.10 kU/L	No significant level detected	Class 0	Allergen results of 0.10-0.34 kU/L for whole peanut are intended for specialist use as the clinical relevance is undetermined. Even though increasing ranges are reflective of increasing concentrations of allergen-specific IgE, these concentrations may not correlate with the degree of clinical response or skin testing results when challenged with a specific allergen. The correlation of allergy laboratory results with clinical history and in vivo reactivity to specific allergens is essential. A negative test may not rule out clinical allergy or even anaphylaxis.
0.10 - 0.34 kU/L	Clinical relevance undetermined	Class 0/1	
0.35 - 0.70 kU/L	Low	Class 1	
0.71 - 3.50 kU/L	Moderate	Class 2	
3.51 - 17.50 kU/L	High	Class 3	
17.51 - 50.00 kU/L	Very high	Class 4	
50.01 - 100.00 kU/L	Very high	Class 5	
Greater than 100.00 kU/L	Very high	Class 6	

**Peanut Allergen Components**



Patient: [REDACTED]  
 ARUP Accession: 22-262-109871

# Allergen, Food, Peanut Reflex Components IgE (Reflex for 3002253 PEANUT R-Not orderable by clients)

Patient: [REDACTED] | Date of Birth: [REDACTED] | Sex: F | Physician: TEST  
Patient Identifiers: [REDACTED] | Visit Number (FIN): [REDACTED]

## Interpretive Information

Severe Risk: Ara h 1, 2, 3, 6, and 9 are the peanut proteins most linked to systemic allergic reactions. If any one of these five components is positive, the patient has a high risk of systemic allergic reaction when exposed to peanuts, regardless of the whole peanut allergen result.

