

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 10/14/1998  
**Gender:** Male  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 00/00/0000 00:00

**Brivaracetam Quantitative, Serum or Plasma**

ARUP test code 3004978

Brivaracetam Quantitative, Serum/Plasma

None Det mcg/mL  
Serum or Plasma  
Reporting Limit: 0.10 mcg/mL

Synonym(s): Briviact(R)  
The recommended steady-state brivaracetam plasma concentration for seizure control is 0.2 to 2.0 mcg/mL. Analysis by High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)  
This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the US Food and Drug Administration.  
Testing performed at NMS Labs, Inc.  
200 Welsh Road  
Horsham, PA 19044-2208  
CLIA 39D0197898

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Brivaracetam Quantitative, Serum/Plasma	22-122-103651	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

**H=High, L=Low, \*=Abnormal, C=Critical**

Unless otherwise indicated, testing performed at: