

HOTLINE: Effective August 19, 2019

0013005

Antibody ID RBC Prenatal-Reflex to Titer

IRL-ABID

Specimen Required: Collect: 3 (7mL) Lavender (K₂EDTA) or Pink (K₂EDTA) AND 1 (7mL) Plain Red.

Specimen Preparation: Do not freeze. Transport 3 (7 mL) whole blood EDTA AND 1 (7 mL) whole blood plain red. (Min: 10 mL EDTA and 3 mL Plain Red)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Separator tubes.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 1 week; Frozen: Unacceptable

Note: This test is for prenatal patients only. Includes: ABO/Rh Type, Rh Phenotype, Direct Coombs, RBC Antibody Identification by various methods. Titers will be performed, at an additional charge, on prenatal specimens for clinically significant antibodies.