

Quarterly HOTLINE: Effective **August 20, 2018**

0050379	Herpes Simplex Virus Type 1 Glycoprotein G-Specific Antibody, IgG by ELISA, CSF	HERPICSF
----------------	--	-----------------

Specimen Required: Collect: CSF.

Specimen Preparation: Transfer 0.5 mL CSF to an ARUP Standard Transport Tube. (Min: 0.2 mL)

Storage/Transport Temperature: Refrigerated.

Remarks: Indicate source on test request form.

Unacceptable Conditions: Specimen types other than CSF. Contaminated, heat-inactivated, **or hemolyzed** specimens.

Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 2 weeks; Frozen: 1 year