

HOTLINE: Effective November 14, 2022

## 0051368 RhD Gene (*RHD*) Copy Number

RHD

Specimen Required: Collect: Fetal genotyping: Amniotic fluid

**OR\_Cultured amniocytes**: Two T-25 flasks at 80 percent confluency.

OR cultured CVS: Two T-25 flasks at 80 percent confluency.

If the client is unable to culture, order test Cytogenetics Grow and Send (ARUP test code 0040182) in addition to this test and ARUP will culture upon receipt (culturing fees will apply). If you have any questions, contact ARUP's Genetics Processing at 800-522-2787 ext. 3301.

WITH maternal cell contamination specimen (see Remarks): Lavender (EDTA), pink (K<sub>2</sub>EDTA), or yellow (ACD Solution A or B).

Parental genotyping: Lavender (EDTA), pink (K<sub>2</sub>EDTA), or yellow (ACD Solution A or B).

Specimen Preparation: Amniotic fluid: Transport 10 mL amniotic fluid in a sterile container. (Min: 5 mL)

Cultured amniocytes AND cultured CVS: Transport two T-25 flasks at 80 percent confluency filled with culture media. Backup cultures must be retained at the client's institution until testing is complete.

Maternal cell contamination specimen: Transport 3 mL whole blood (Min: 1 mL)

Whole blood (parental genotyping): Transport 3 mL whole blood (Min: 1 mL)

<u>Storage/Transport Temperature:</u> Amniotic fluid, cultured amniocytes and cultured CVS: CRITICAL ROOM TEMPERATURE. Must be received within 48 hours of shipment due to lability of cells.

viust be received within 48 nours of snipment due to lability of cells.

Whole blood or maternal cell contamination specimen: Refrigerated.

<u>Remarks:</u> Maternal specimen is recommended for proper test interpretation if contamination of the fetal specimen from the mother is suspected. Order Maternal Cell Contamination. Patient History Form is available on the ARUP website or by contacting ARUP Client Services.

Unacceptable Conditions: Frozen specimens in glass collection tubes.

Stability (collection to initiation of testing): Amniotic fluid, cultured amniocytes and cultured CVS Fetal Specimen: Ambient: 48 hours; Refrigerated: Unacceptable; Frozen: Unacceptable

Whole blood or maternal cell contamination specimen: Ambient: 72 hours; Refrigerated: 1 week; Frozen: 1 month

## Interpretive Data: Refer to report

Kelei to lepon

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA-certified laboratory and is intended for clinical purposes.

Counseling and informed consent are recommended for genetic testing. Consent forms are available online.