

TEST CHANGE

Phenol Exposure Quantitative, Urine

0091260, PHENOL U

Specimen Requirements:

Patient Preparation:

Collect: Urine.

Specimen Preparation: Transfer 4 mL urine to an ARUP **standard transport tube**~~Standard Transport Tube~~. (Min: 1.9 mL) Preservative-free urine specimens are recommended. Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions: Specimens preserved with **benzoic acid**~~Benzoic Acid~~.
Specimens received room temperature.

Remarks:

Stability: Ambient: 4 days; Refrigerated: 1 week; Frozen: 22 months

Methodology: Quantitative Gas Chromatography/Colorimetry

Performed: Varies

Reported: **8-16**~~5-12~~ days

Note:

CPT Codes: 84600; 82570

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report