

## TEST CHANGE

### Phenol Exposure Quantitative, Urine

0091260, PHENOL U

#### Specimen Requirements:

##### Patient Preparation:

Collect: Urine.

Specimen Preparation: Transfer **24** mL urine to an ARUP standard transport tube. (Min: **1.9** mL) Preservative-free urine specimens are recommended. Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: **Frozen**~~Refrigerated~~. Also acceptable: ~~Refrigerated~~**Frozen**.

Unacceptable Conditions: Specimens preserved with benzoic acid. Specimens received room temperature.

##### Remarks:

Stability: Ambient: 4 days; Refrigerated: 1 week; Frozen: 22 months

Methodology: Quantitative Gas Chromatography-~~Mass Spectrometry (GC-MS)~~  
/ Colorimetry

##### Note:

CPT Codes: 84600; 82570

New York DOH Approval Status: This test is New York DOH approved.

##### Interpretive Data:

##### Reference Interval:

~~Refer to~~**By** report