
0099169 JC Virus by PCR

JC VIRUS

Specimen Required: Collect: Lavender (EDTA), pink (K₂EDTA) or serum separator tube. OR CSF or urine.

Specimen Preparation: Separate serum or plasma from cells. Transfer 1 mL serum, plasma, CSF or urine to a sterile container. (Min: 0.5 mL)

Storage/Transport Temperature: Frozen.

Remarks: Specimen source required.

Unacceptable Conditions: Heparinized specimens, **tissues in optimal cutting temperature compound.**

Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 5 days; Frozen: 30 days