

---

**2000270 Isohemagglutinin Titer, IgM**

**IRL ISO M**

**Specimen Required:** Collect: Lavender (K<sub>2</sub>EDTA), or Pink (K<sub>2</sub>EDTA).

Specimen Preparation: **Do not freeze red cells.** Transport 7 mL whole blood. (Min: 3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Separator or gel tubes.

Stability (collection to initiation of testing): Ambient: **Unacceptable**; Refrigerated: 1 week; Frozen: Unacceptable