
2000280 Isohemagglutinin Titer, IgG and IgM

IRL ISO MG

Specimen Required: Collect: Lavender (K₂EDTA) or Pink (K₂EDTA).

Specimen Preparation: **Do not freeze.** Transport 14 mL whole blood. (Min: 6 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: **Separator or gel tubes.**

Stability (collection to initiation of testing): Ambient: **Unacceptable**; Refrigerated: 1 week; Frozen: Unacceptable