

HOTLINE: Effective February 16, 2021

**2002109 Protein Electrophoresis with Reflex to Immunofixation, Serum SPEP REFLEX**

**Methodology:** Quantitative Capillary Electrophoresis/Qualitative Immunofixation Electrophoresis/Quantitative Immunoturbidimetry/Quantitative Spectrophotometry

**Reference Interval:**

Test Number	Components	Reference Interval		
0050640	Protein Electrophoresis, Serum	Effective August 19, 2019		
		Test Number	Components	Reference Interval
			Total Protein, Serum	Refer to report
			Albumin	Refer to report
			Alpha-1 Globulins	Refer to report
			Alpha-2 Globulins	Refer to report
			Beta Globulins	Refer to report
			Gamma	Refer to report
0050350	Immunoglobulin G	Effective February 16, 2021		
		0-2 years: 242-1108 mg/dL 3-4 years: 485-1160 mg/dL 5-9 years: 514-1672 mg/dL 10-14 years: 581-1652 mg/dL 15-18 years: 479-1433 mg/dL 19 years and older: 768-1632 mg/dL		
0050340	Immunoglobulin A	Effective February 16, 2021		
		0-2 years: 2-126 mg/dL 3-4 years: 14-212 mg/dL 5-9 years: 52-226 mg/dL 10-14 years: 42-345 mg/dL 15-18 years: 60-349 mg/dL 19 years and older: 68-408 mg/dL		
0050355	Immunoglobulin M	Effective February 16, 2021		
		0-2 years: 21-215 mg/dL 3-4 years: 26-155 mg/dL 5-9 years: 26-188 mg/dL 10-14 years: 47-252 mg/dL 15-18 years: 26-232 mg/dL 19 years and older: 35-263 mg/dL		