
2012278

Carisoprodol Urine with Reflex to Quantitation

CARIS RFXU

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Breast milk. Pharmaceutical preparation. Specimens exposed to multiple freeze thaw cycles.

Stability (collection to initiation of testing): Ambient: **1 week**; Refrigerated: **1 month**; Frozen: **1 month**