

HOTLINE: Effective February 22, 2022

3003751 JAK2 (V617F) Mutation by ddPCR, Quantitative

JAK2V617FQ

Specimen Required: Collect: Whole blood or bone marrow: Lavender (EDTA).

Specimen Preparation: Whole Blood: Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)

Bone Marrow: Do not freeze. Transport 3 mL bone marrow. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or fresh or frozen tissue. Specimens collected in anticoagulants

other than EDTA or sodium heparin. Clotted or grossly hemolyzed specimens.

Stability (collection to initiation of testing): Refrigerated: 7 days; Frozen: Unacceptable