

HOTLINE: Effective August 15, 2022

---

3004383

**Cerebral Autosomal Dominant Arteriopathy with Subcortical Infarcts and Leukoencephalopathy, CADASIL (*NOTCH3*), Sequencing**

**NOTCH3 NGS**

**Specimen Required:** Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B).

Specimen Preparation: Transport 5 mL whole blood. (Min: 3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue.

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable