

HOTLINE: Effective July 5, 2022

3004547 Beta Globin (*HBB*) Sequencing

BG NGS

Specimen Required: Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B).

New York State Clients: Lavender (EDTA)

Specimen Preparation: Transport 3 mL whole blood. (Min: 2 mL)
New York State Clients: Transport 7 mL whole blood (Min. 3 mL)

Storage/Transport Temperature: Refrigerated

Remarks: Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue. Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

Private Information Page 1