
3004716

Galactosemia (*GALT*) Sequencing and Deletion/Duplication

GALT NGS

Specimen Required: Collect: Lavender or pink (EDTA) or yellow (ACD solution A or B).

New York State Clients: Lavender (EDTA).

Specimen Preparation: Transport 3 mL whole blood. (Min: 2 mL)

New York State Clients: Transport 3 mL whole blood. (Min. 3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue.

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable